**DAYVIA**<sup>®</sup>  
LUMINOTHERAPIE

---

Lunettes de luminothérapie

**SunActiv/02**

Lunettes de luminothérapie

## SunActiv/02

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P12: User guide (EN)

P21: Bedienungsanleitung (DE)

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## Lunettes de luminothérapie

# SunActiv/02

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
## Light therapy glasses

# SunActiv/02

Figure 1: The SUNACTIV telescope and its accessories

## INTRODUCTION

Congratulations for buying your DAYVIA light therapy glasses.

Used well, these glasses will help you regain the light you are missing.

The lack of light shows itself in some people by several symptoms:

-Irritable, and bad humor.

-Significant increase in appetite.

-Reduction in activity, and productivity in work linked to lack of energy.

-Trouble sleeping, and chronic fatigue.

-Sexuality disorders (decreased libido)

Numerous international studies show that light therapy helps to regain energy, shape, ability to concentrate and a good mood such as when you are sunbathing in the summer.

The light therapy glasses also helps to adjust the sleeping / waking rhythm for people who have staggered work schedules (night workers), people who suffer from phase delay (late falling asleep and difficulty waking up in the morning), or advanced phase (early wake up) or people who travel and cross several time zones (jet-lag).

The effect of light depends on its intensity but especially on its spectral curve (color). Indeed, certain receptors located in our eyes are very sensitive to light with certain wavelengths around 500 nanometers (blue, cyan and / or green), they have a major role on our biological clock.

Thus, a cyan monochromatic light included in this spectrum such as that emitted by this medical device can prove to be as effective as a much more intense white light.

Please read these instructions carefully to take advantage of the advantages of your DAYVIA light therapy glasses.

## SECURITY MEASURES

- The device must be installed and used in accordance with the information provided in the accompanying documents.

- Do not use the glasses in a dangerous situation such as driving a vehicle.

- Do not use the DAYVIA glasses in a humid environment (eg in a bathroom, near a bathtub or shower, outside ...).

- Keep the charger's power cord away from hot spots.

- When used in a warm environment (ambient temperature of 40 ° C), certain parts of the glasses body can reach a T ° of 43 ° C to the touch, it is necessary to limit the contact of these parts with the skin to 10 minutes in this hot environment. The parts which normally come

Hinto contact with the face (nasal support and arms) are not affected, it is therefore sufficient to position the glasses well so that it can be used without any time limit in a hot environment.

- Use the charger supplied with your light therapy product or any system compatible with the characteristics specified in these instructions, the assembly forms an electro-medical device.

- The mains plug serves as a power disconnect switch. The appliance must be recharged with an easily accessible mains socket.

- These glasses should not be used by children or persons with reduced physical or intellectual capacities.

- Any modification of the device is forbidden.

- The charger must not be used if the power cord and / or the power supply are apparently damaged or defective.

- In accordance with the regulations for medical devices:

- The device needs special precautions, it can be used in all premises including domestic premises and away from RF sources where the intensity of electromagnetic disturbances is high.

- Avoid using this unit next to other equipment because this may cause malfunction. If this is necessary, this appliance and other appliances should be observed for normal operation.

- Do not use RF portable communication devices (including peripherals such as antenna cables and external antennas) closer than 30 cm (12 inches) to any part of the DAYVIA glasses. Otherwise, the performance of these devices may be impaired.

## PRECAUTIONS

The glasses must be used in an environment illuminated by another light source to ensure optimum comfort in the eyes and to avoid too much contrast between the glasses and the eyes. The glasses are complementary to normal lighting (electric or natural light).

Unlike daylight or other lamps on the market, this medical light therapy device emits no harmful radiation such as UV for example.

Studies have shown that light therapy can be practiced safely.

However, before starting your treatment, please consult your doctor if:

- You suffer or have suffered from depression.

- You are hyper sensitive to light.

- You have an eye disease.

- You suffer from mental disorders.

- You are receiving treatment with psychotropic or

photosensitizing medications.

If you have any problems (such as headaches or eye strain) as a result of using the light therapy glasses please contact your doctor. It is often enough to begin treatment with shorter sessions (10 minutes) and gradually increase them each day to the recommended duration.

## COMMISSIONING / RECHARGING / SHUTDOWN

The DAYVIA glasses are carefully checked before being packaged and shipped.

You will find in your box: A user manual, a charger with its cord and your light therapy glasses in their case.

### **Before your first use we recommend that you charge your glasses:**

1. Connect the standard USB cable connector to the supplied charger port, and then insert the Micro USB connector into the socket on the glasses (see Figure 1);
2. Plug the charger into a wall outlet;
3. When the bezel is fully charged, unplug the bezel from the charger.

Information: To optimize the life of the battery we recommend to carry out a complete maintenance charge every 3 months when not in use.

### **Charge indication of your glasses:**

Your glasses have 6 LEDs which light up gradually in groups of 2 at low intensity and continuously during charging.

- When you start the charge, you will see lighting first of the 2 LEDs on the far left, then the 2 middle LEDs and finally the 2 LEDs on the far right, each time symbolizing the charging of the third of the total load

- When the charge is complete, these 6 LEDs remain on continuously.

From this moment, you can unplug your glasses from the charger and use it for your session.

- In case of internal problem, the 2 LEDs on the right flash very quickly at high intensity for 5 seconds. .

- Your glasses are not operational during charging.

Information: The charging of your glasses implies a normal heating (about 50 ° C for the battery at the end of charge).

### **Autonomy of your glasses:**

The autonomy of the well-charged device varies between 3 and 5 hours, depending on the level of intensity chosen during your sessions.

When starting the glasses, after the 3 intensity levels, the 4th press on / off allows you to view the charge level and make sure that it is sufficient to start a session.

When you turn off your bezel you will also be able to see if there is sufficient autonomy for your next session.

After your pressing for 1 second on the OFF button of the device, the bezel switches off and:

- Either the 6 LEDs reappear and light up at low intensity for 3 seconds = the remaining autonomy is greater than 90%, therefore sufficient for your next session.

- Either 4 LEDs reappear and light up at low intensity for 3 seconds = the remaining battery life is over 40%, so enough to carry out your session the next day.

- Either 2 LEDs reappear and light up at low intensity for 3 seconds = the remaining battery life is between 10 and 40% so you have to recharge it before using it again.

- - If all the LEDs suddenly turn off during your session = the remaining battery life is less than 10%. In this case, the 2 LEDs on the left flash at low intensity for 3 seconds and turn off automatically. You need to charge your glasses.

### **Switching on and off:**

Make sure your glasses are fully charged to complete the entire session.

Pressing the ON / OFF button briefly switches the bezel on at intensity level 3 (maximum intensity).

A second short press causes a decrease in the intensity at level 2.

A third press triggers level 1 (minimum intensity).

A fourth press indicates the charge level

Once Level 1 is reached, you can increase the intensity by pressing again and the cycle begins again : 1 press to return to level 3, then another to go to reach level 2, then level 1 and after the load level. Further short presses activate the cycle again and so on.


To turn off your glasses, press and hold the ON / OFF button (1 second). The charge indication informs you if the autonomy is sufficient to carry out your session of the following day (see paragraph: Autonomy of your glasses).

## **RECOMMENDATIONS FOR USE**


### **Positioning of the glasses**

Figure 2: How to position my glasses correctly

- If you do not wear glasses

- If you wear glasses:

The ideal position is that which allows you to receive a maximum of light at the level of the eyes (see figure 2), **we advise you to observe yourself in front of a mirror when you use it for the first time to make a good adjustment of the nose rest.**

The nose rest is adjustable with 3 notches but it can also be removed to reverse its position and thus optimize the diffusion according to each face morphology and whether or not you wear spectacles.

If you do not wear glasses, the ideal position is in most cases is to opt for the nose rest facing down.

If you are wearing eyeglasses, the nosepiece of your DAYVIA bezel should be above the frame of your spectacles and for the majority of the morphologies it is recommended that the nose rest be pointing upwards and not downwards.

Do not try to fix with your eyes the emission points of the light above the eyes, there is no danger. You can continue your activities throughout your session (breakfast, reading, office work, computer, walking, exercise bike ..).

For a pleasant, glare-free use, always make sure that the room is lit (especially for the first days of use) to minimize the contrast between the illumination of the light therapy device and its environment.

### **Recommendations on frequency, timing and duration of use:**

It is advisable to start the light therapy sessions as soon as you feel a lack of light before the first symptoms of the winter blues appear.

In prevention or treatment of winter blues, it is preferable to practice the sessions in the morning after waking up.

Make a session each morning in a 3-week cycle during the dark period. These figures must be adapted to each person according to the results obtained after a few days of treatment.

The table below shows the recommended daily duration according to the intensity level chosen:

<b>Level (at start-up)</b>	<b>Recommended time (in minutes)</b>
Level 3	30
Level 2	45
Level 1	60

Once the session is over, the 6 LEDs flash for 3 seconds and turn off automatically.

If you have the opportunity to practice at least 15 minutes outside during the day, even in winter, you can reduce your daily light therapy session to 20 min at intensity level 3 or 30 min at level 2 and 45 min at level 1.

In general, all people wishing to compensate for their lack of light (winter season, confined work ...) by light therapy should favor an exposure in the morning. This morning exposure, ideally just after waking up, is also recommended for people who are late in life (late falling asleep and difficulty emerging in the morning) and people preparing 2 to 3 days before traveling east (crossing more than 3 time zones), the effects of the jet lag will be minimized.

Light therapy at the end of the day is reserved for people working at night, for people who are early in life (sleeping too early and getting up too early) and for people

preparing for a trip to the west (crossing of more than 3 time zones) to avoid jet lag.

Your glasses are also suitable when you have a slight fatigue at a time in the day, for example after lunch. The light will thus slow down or block the secretion of melatonin whose secretion reduces your alertness.

## MAINTENANCE OF THE APPLIANCE

Your glasses are recommended for personal use only and we recommend that you place it in their case in a dry place after each use. Never use abrasives or detergents to clean the appliance: Use a slightly moistened cloth, taking care that no liquid enters the appliance.

The battery is not replaceable. Its life span, like that of the light therapy glasses, is at least 5 years under normal use, there is no maintenance. The life of the charger supplied with the glasses is at least 3 years.

## CHARACTERISTICS OF THE DEVICE

Specifications are subject to change without notice.

- Internal power supply of the glasses: Rechargeable battery AAA NiMH 1.2 V 1000 mAh
- Operating temperature: between + 5 ° C and + 40 ° C
- Humidity rate during use: between 15% and 93% (without condensation)
- Penetration of foreign bodies and water for the glasses: IP22
- Atmospheric pressure during use: between 700 hPa and 1060 hPa
- Storage and transport temperature: between -25 ° C and + 70 ° C with relative humidity up to 90% without condensation and water vapor pressure not exceeding 50 hPa
- Overvoltage category II
- Size of the glasses: 14.1 x 1.9 x 18.2 (W x H x D in cm)
- Weight of the glasses without charger: 50 g
- Weight of the glasses + charger and cord: 80 g
- Class IIa medical device (CE0459)

DAYVIA Sun Activ has been certified to comply with all standards for medical devices in Europe:

- conforms to electrical security.
- conforms to photobiological risks of light : group 0 (no risk).
- conforms to EMC (Electromagnetic compatibility): Class B and Group 1 (Emission RF CISPR 11). Electrostatic discharge immunity test: 15kV for air discharge, 8 Kv for contact discharge. Radiated, radio-frequency,

electromagnetic field immunity test: 10V/m. Electrical fast transient/burst immunity test: 2 kV. Surge immunity test:  $\pm$  1 kV differential mode and  $\pm$  2 kV common mode. Immunity to conducted disturbances, induced by radio-frequency fields: 150kHz to 80MHz 3V, 6V for other frequency bands. Power frequency magnetic field immunity test: 3A/m. Voltage dips, short interruptions and voltage variations immunity tests: 0% 10ms, 70% 500ms, 0% 5000ms.

### Charger specifications supplied with the bezel:

- Manufacturer: Shen Zhen Bizoe Electronic Technology
- Model: BZ-A011
- Charger power supply: 0.15 A
- Charger voltage: 100-240 V AC
- Output voltage of the charger: 5 V
- Frequency of charger: 50-60 Hz
- Insulation of the charger: Class II
- USB cable length of charger: 1m00

## GUARANTEE AND SUPPORT

If you need information or support, please visit the DAYVIA website at <http://www.dayvia.com/>.

## EXPLICATION DES SYMBOLES

Consult the supplied user's guide.

Manufacturer

Serial identification.

This appliance includes applied parts of type BF (left and right branches + nose rest of the bezel)

This device must not be discarded with house-hold wastes, take it to a waste disposal facility where it can be recycled. By doing so, you will contribute to the protection of the environment.

IP22 Degree of protection against liquids and particles.

Temperature limits .

Humidity limits.

This device is certified to conform to the guidelines regulating medical devices.

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