lumie® Bodyclock™

ACTIVE 250

USER GUIDE

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BRUKSANVISNING

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Safety

If the product is damaged please contact us before use.

Do not use Bodyclock Active if the lamp cover is broken, damaged or missing. If you need to remove the lamp cover, unplug the unit first and let it cool down.

The lamp cover gets warm when the light is on, so keep the unit clear of anything that may be damaged by heat.

Do not open the base – there are no user-serviceable parts inside.

For indoor use only.

Keep away from water and damp.

Keep the cable away from heated surfaces. If the cable is damaged it must be replaced by Lumie, a Lumie distributor or suitably qualified person in order to avoid a hazard.

Not to be used by or on individuals who are unconscious, anaesthetised or incapacitated e.g. under the influence of alcohol.

Electromagnetic compatibility

This product conforms to the electromagnetic compatibility standard for medical devices and is unlikely to cause interference affecting other electrical and electronic products in its vicinity. However, radio transmitting equipment, mobile phones and other electrical/electronic equipment can interfere with this product; if this happens simply move the product away from the source of interference. In rare cases where product function is lost due to EMC effects, disconnect the product from power source for 5 seconds and then restart. If the problem persists, please contact support.

For further detailed information on the electromagnetic environment in which the device can be operated refer to the website: www.lumie.com/pages/e-m-c.

Disposal

At the end of its life, do not throw the unit out with normal domestic waste. Please take it to an official recycling point to be disposed of responsibly.

 SN	Serial number		Do not place in household waste. Dispose of at a recycling centre.		Fuse
	European Declaration of Conformity: Class I Medical Device		Read instructions for use		Double-insulated (Class II electrical)
	Manufacturer: Lumie, 3 The Links, Bar Hill, Cambridge, CB23 8UD, UK.		Caution		AC (alternating current) power supply
	Keep dry		Caution – hot surface		Member of the 'Green Dot' recycling scheme
	For indoor use only				Recyclable

Introduction

Lumie Bodyclock Active wakes you with a sunrise that naturally helps to regulate your sleep and wake patterns. You'll wake up in a better mood, find it easier to get up, and feel more energetic during the day. There's also a fading sunset to help you unwind at the end of the day and drift off to sleep.

Ideally, you should use Bodyclock Active all year round to help keep your sleep and wake cycles on track.

In clinical trials, 'dawn simulators' such as Bodyclock Active have been shown to improve mood, energy, productivity and the quality of sleep and awakening[†]. They can also help sufferers of SAD (Seasonal Affective Disorder) and winter blues cope with dark winter mornings.

Lumie was the first to launch a dawn simulator range and began selling an early version of Bodyclock in 1993.

Please read these instructions carefully to get the most out of your Bodyclock Active.

[†] Acta Psychiatr Scand 1993 Jul; 88 (1):67-71.

Features

Bodyclock Active is supplied ready to use with the settings we've found most popular with users; all you need to do is fit the bulb, set the current time and set the time you want to wake up.

The standard settings are:

- An optional 30-minute sunset to help you get to sleep.
- A 30-minute sunrise to wake you gently, followed by an alarm beep.
- A clock display that automatically dims to a low level during sunset.

Alternatively, Bodyclock Active provides a range of other options to suit your personal preferences. For example, you can choose:

- A longer or shorter sunrise and sunset.
- A different alarm sound or the radio.
- A clock display that is off or stays bright at night.

Getting started


Here's how to set up your Bodyclock Active for the first time.

The bulb is supplied packed in a separate box and needs to be fitted into the Bodyclock Active before use.

The lamp cover may be supplied already separated from the base, in which case skip the next step.

To remove the lamp cover

The lamp cover is attached to the base of the Bodyclock Active by two internal clips, one at either end of the cover:

- Support the base of the Bodyclock Active with one hand.
- Push firmly at one end of the lamp cover to release the clip.
- Remove the lamp cover.

To fit the bulb


- Screw the bulb into the holder.
- Replace the lamp cover by pressing it down until the clips at each end pop into place.

To position Bodyclock Active

- Decide where to put your Bodyclock Active.

The ideal place is on a bedside table near the head of your bed, near a mains socket.

To set up Bodyclock Active

- Plug in the unit and switch on at the mains.

You will then be prompted to enter the essential settings:

Display shows (for example)	What you do
SET TIME 13:00	Use \ominus or \oplus to set the current time (24-hour format). Press \odot to continue.
SET WAKEUP TIME 7:00	Use \ominus or \oplus to set the time you want to be woken up (24-hour format). Press \odot to continue.
WAKE SOUND 15:17	At this stage you can ignore the remaining options. Wait 12 seconds for the normal display.

To activate the alarm

When you're ready to go to sleep:

- Press \odot .

The display shows \odot to confirm that the alarm is set with a sunrise and shows the alarm type and wake-up time you've set; for example:

ALARM BEEP
 \odot 8:15

Using the bedside light and radio

Your Bodyclock Active can be used as a bedside light which you can set to your preferred brightness.


Bodyclock Active also incorporates an FM radio and has five other sound options to help you relax or mask other sounds when you're trying to sleep:

Dawn chorus	Rooster	Waves	White noise	Alarm beep	Radio
					

You can also use the radio or any of the sounds as an alarm sound.

Note: Extend the wire aerial for best reception.


To use the bedside light





- Hold down ⊕ to set the light level you want.
- Hold down ⊖ to turn the light off.

Note: If you leave the light on it will automatically turn off after 90 minutes.

To use the radio





- Press  to turn the radio on.
- Use  and  to tune to a station.
- Press  or  to adjust the volume.

A short press tunes manually. A longer press skips to the next station when you release the button.

- Press  to turn the radio off.

Note: If the radio has been on for less than 12 seconds press  again to skip the sounds.


To play a well-being sound

- Press  twice to select the sounds.
- Use  and  to select the sound you want.
- Press  or  to adjust the volume of the sounds.
- Press  to turn the sounds off.

Using the alarm

When you're ready to go to sleep you need to explicitly activate the alarm for the following morning. This starts dimming the light to give a sunset.

Bodyclock Active will wake you up with a gentle sunrise, ending with an alarm sound at the time you have set.

To activate the alarm

When you're ready to go to sleep:

- Press .

The display shows  to confirm that the alarm is set with a sunrise, and shows the alarm sound and wake-up time that you've set; for example:

ALARM BEEP
 8:15 

 is shown if you have chosen one of the VOLUME FADE AUTO or VOLUME FADE PERM settings to fade off the radio or well-being sound; see *Choosing sleep options*, page 12.

If you leave the bedside light on, it will fade off gradually to help you go to sleep. Alternatively, to turn off the light immediately:

- Hold down .

To snooze

- Press .


This gives you another nine minutes before the sunrise (and alarm sound, if you've chosen one) is repeated.

To cancel the alarm sound

- Press .

To turn off the light

When the alarm is off:

- Hold down .

Alternatively, the alarm sound and light will turn off automatically after 30 minutes.

Choosing wake-up options

You can choose whether you wake to a sunrise followed by one of the relaxation/well-being sounds, to a sunrise and the radio, to a simple sunrise only or to a sunrise ending with a flashing light.

To set a wake-up option

- Press and hold  until SET WAKEUP TIME shows.
- Press  twice until the current WAKE setting is shown; for example: WAKE SOUND.

To choose a sound alarm

- Press  or  to set the WAKE setting to WAKE SOUND.
- Press .
- Press  or  to step through and select the sound option you want:

Dawn chorus	Rooster	Waves	White noise	Alarm beep
				


- Press  to set the sound option.

You are then prompted to set the wake-up volume.

- Press  or  to change the volume.
- To set the volume press  or wait 12 seconds.

To choose a radio alarm

- Press  or  to set the WAKE setting to WAKE RADIO.
- Press .
- Choose an FM radio station by pressing  (to tune down) or  (to tune up).

Note: Don't use the usual radio tuning buttons  or .

- Press  to set the station.

- Press  or  to change the volume.
- To set the volume press  or wait 12 seconds.

To choose a sunrise-only alarm

- Press  or  to set the WAKE setting to WAKE LIGHT.
- Press  or wait 12 seconds.

This wakes you with a simple sunrise without any sound.

To choose a flashing light alarm

- Press  or  to set the WAKE setting to WAKE FLASH.
- Press  or wait 12 seconds.

This gives you a sunrise followed by the light flashing on and off at the wake-up time without any sound.

Choosing sleep options

You can choose whether the time display is bright, dimly illuminated or off at night to minimise sleep disturbance.

If you go to sleep listening to the radio or a well-being sound you can choose whether to have the volume fade down gradually to silence or to a low level until you wake up (for example, to block out unwanted noise).

To set the display brightness

- Press and hold  until SET WAKEUP TIME shows.
- Press  repeatedly until the current DISPLAY option is shown; for example, DISPLAY AUTO.
- Press  or  to step through and select a DISPLAY setting:

Option	What it means
DISPLAY AUTO	The display dims to a low level when the light is low or off.
DISPLAY HIGH	Keeps the display at its brightest setting without dimming at all.
DISPLAY AUTO OFF	The display dims until it is completely off when the light is low or off.

- To set, press  or wait 12 seconds.

Note: If you have set DISPLAY AUTO OFF and you need to see the time when the display is off, press  to light the display. After 5 seconds it will revert to being off.

To fade the radio or sounds

- Press and hold  until SET WAKEUP TIME shows.
- Press  repeatedly until the current VOLUME FADE option is shown; for example, VOLUME FADE OFF.
- Press  or  to step through and select a volume-fade option:

Option	What it means
VOLUME FADE OFF	The volume doesn't fade.
VOLUME FADE AUTO	The volume fades to silence over the duration of the sunset.
VOLUME FADE PERM	The volume fades over the duration of your sunset and then stays on at a low level.

- To set, press  or wait 12 seconds.

Note: If you have set VOLUME FADE AUTO or VOLUME FADE OFF,  is shown in the display to remind you.

Customisation options

This section summarises the wake-up and sleep options already described, and explains the other options available in Bodyclock Active.

To check a setting

If the alarm is activated (☀️ shown on display):

- Press  to cancel it.
- Hold down  until the display shows the first setting:

SET WAKEUP TIME

- Press  to step through the options until the one you want is shown; see table below.


To change a setting

When the setting you want is shown:

- Press  or  to step through the alternative options.
- Press  to confirm your setting.
- Either press  to choose another setting or wait 12 seconds to return to normal operation.

The following table explains the meaning of each setting:

Default display	What it does	Options
SET WAKEUP TIME	Sets the time you want to be woken up.	Example: 07:00.
SET TIME	Sets the current time.	Example: 13:00.
WAKE SOUND	Sets the type of alarm; see <i>Choosing wake-up options</i> , page 10.	WAKE SOUND WAKE RADIO WAKE LIGHT WAKE FLASH

Default display	What it does	Options
ALARM BEEP	When WAKE SOUND is selected, sets the sound to use as the alarm sound; see <i>Choosing wake-up options</i> , page 10.	ALARM BEEP DAWN CHORUS ROOSTER WAVES WHITE NOISE
	When WAKE SOUND or WAKE RADIO is set, allows you to set the wake-up volume for your alarm sound or radio using the  or  buttons.	
DISPLAY AUTO	Sets whether the display dims; see <i>Choosing sleep options</i> , page 12.	DISPLAY AUTO DISPLAY HIGH DISPLAY AUTO OFF
SUNRISE 30	Sets the sunrise duration, in minutes.	SUNRISE 30 45 60 90 15 20
SUNSET 30	Sets the sunset duration, in minutes.	SUNSET 30 45 60 90 15 20
VOLUME FADE OFF	Sets whether the radio or sounds fade off during sunset; see <i>Choosing sleep options</i> , page 12.	VOLUME FADE OFF VOLUME FADE AUTO VOLUME FADE PERM
NIGHTLIGHT OFF	Sets an optional nightlight, which stays on at a low level all night.	NIGHTLIGHT OFF NIGHTLIGHT ON
SECURITY OFF	Sets security mode; see <i>To set security mode</i> , page 16.	SECURITY OFF SECURITY ON
ENGLISH	Sets the language for the display messages.	ENGLISH FRANÇAIS DEUTSCH NEDERLANDS SVENSKA SUOMI ESPAÑOL

To have a sunset without activating the alarm

If you don't need to be woken by the alarm, for example if the next day is a weekend, you can still activate the sunset feature to help you go to sleep:

- Press both  and  at the same time.

The display shows SUNSET and the light slowly starts to fade.

This also dims the display and fades the radio/well-being sound if you've selected these options.

To cancel the sunset and return to a normal display:

- Press both  and  at the same time again.

To set security mode

Security mode allows you to use Bodyclock Active as a security light when you're away. The light switches on and off at random times between 16:00 and 23:00 to give the impression that there's someone at home.

- Press and hold  until SET WAKEUP TIME shows, then press  repeatedly until SECURITY OFF shows.
- Press  to select SECURITY ON.
- Wait 12 seconds.

In security mode you cannot use any of the other functions, and between the hours of 16:00 and 23:00 the display shows:

SECURITY

To cancel security mode and return to normal operation:

- Press any button.

Frequently asked questions

I tried out the sunrise but the light comes up really quickly

You probably didn't allow enough time for a full sunrise. For example, the default sunrise takes 30 minutes so if the clock shows the time as 11.00 then you need to set the wake-up time to at least 11.30.

What happens when I unplug my Bodyclock Active or there's a power cut?

Your Bodyclock Active keeps time for one hour; other settings are permanent, so you won't need to reprogram the unit, just plug it back in. If you need to adjust the time, see *Customisation options*, page 14.

If you receive a strong FM radio signal, the clock keeps time when tuned to a station with RDS (Radio Data System) as it includes an embedded time signal.

The display is too bright when I'm trying to sleep; how can I change it?

See *To set the display brightness*, page 12.

What can I do if the display is working but the light doesn't come on?

Check to make sure the bulb is firmly screwed into the socket. If this doesn't work, replace the bulb.


Why is nothing happening when I press the buttons?

You may have accidentally set the unit to demonstration mode (only useful for shops, exhibitions etc!) so try pressing  for 12 seconds. If that doesn't work or if the display shows RESET, see *How do I restore the default factory settings?* below.

How do I restore the default factory settings?

Turn off the power at the mains. Press and hold  as you turn the power back on. Release the button when the software version number shows.

What can I do if neither the display nor light are working?

A fuse may have blown. With the unit unplugged, try replacing the 3 ampere  fuse in the plug or the 2 ampere fuse in the unit; simply unscrew the black fuse holder underneath.

Please contact Lumie if you need a replacement bulb or if you still have a problem with your Bodyclock Active.

Maintenance

Cleaning

Clean only with a dry or lightly moistened cloth.

Replacement bulbs

Always use the correct bulb (E14 42 W maximum halogen bulb). Replacements are available from Lumie.

Warranty

This unit is covered under warranty against any manufacturing defect for two years from the date of purchase. This is in addition to your statutory rights. If possible, use the original packaging (including internal pieces) to return your Bodyclock Active for service. Make sure the plug is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first. Please call Lumie for further advice.

Technical description

Subject to change without notice.

Lumie Bodyclock Active 250

230 Vac, 50 Hz, 50 W 

Halogen bulb E14 42 W maximum

Operating temperature range: 5 °C to 35 °C

Storage/transport temperature range: -25 °C to +70 °C

Class I Medical Device 

Contact us

Lumie is Europe's leading specialist in light therapy. Our products are based on the latest scientific studies and we focus solely on light therapy and its benefits.

If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on +44 (0)1920 780500 or email info@lumie.com.

Lumie, 3 The Links, Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, UK.

Lumie is a registered trademark of:

Outside In (Cambridge) Limited, incorporated in England and Wales.

Registered Number: 2647359. VAT No: GB 880 9837 71.

Bodyclock is a registered trademark of Outside In (Cambridge) Limited.

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© Lumie 2015
Bodyclock is a registered trademark of Lumie.
Outside In (Cambridge) Limited met rechtspersoonlijkheid in Engeland en Wales.
Bodyclock is een gedeponeerd handelsmerk van Outside In (Cambridge) Limited.
Outside In (Cambridge) Ltd