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User manual

(climsom)Flex

## Practice of cervical decompression:

1. Place the cervical cushion on a flat and preferably firm surface (carpet for example).
2. Lie down in front of the cervical cushion, first on your elbows.  
You can position yourself on either side of the cushion, depending on what is most comfortable for you:

In this sense, the heat will be felt on the neck and towards the spine.

In this sense, the heat will be felt on the neck and towards the head.

3. Bend your knees and gently position your neck in the curved part of the neck cushion.  
You can extend your arms along your body, or place your hands on each side of your head.  
Flash the QR code to see photos and videos.

4. Stay in this position and concentrate on your breathing.  
Breathe in slowly for 5 seconds, breathe out slowly for 5 seconds (optional: hold your breath for a few seconds after breathing out).

The aim is to let your body adapt naturally to the cervical cushion, a sign of effective cervical decompression.

**Recommended duration of use:** 10 minutes. Adjust this time downwards or upwards according to how you feel.

## Use of heating mode:

The cervical cushion is covered with a heated cover (removable).  
The heating mode will help you release muscle tension and soften your neck, making it easier to decompress your cervical vertebrae.

For more comfort, we advise you to leave the cover on the cushion, even when it is not in use.  
This will avoid contact with the colder material of the cushion.

1. Connect the USB cable from the cover to the supplied 220V EU plug adapter and then to an electrical outlet, or directly to a computer or a power bank.
2. To turn it on, short press the single button on the control unit on the cable.

The cover will turn on by default at the highest temperature (red light : 50°C - 120°F).

Short press to change the temperature :

Green light : 45°C / 115°F - Blue light : 40°C / 105°F

Long press to set an automatic switch-off after 30 minutes.  
Once programmed, the indicator light will flash.

3. To switch off the heating mode, press the button several times until the light is off.  
The number of presses will depend on the temperature being used.

If you experience persistent pain or pain during use, stop and consult a physician or physical therapist. Do not «force», pay attention to your feelings.

**Recommended by François Stévignon, French physiotherapist for more than 30 years, graduate in osteopathy and vertebral therapy, author of the method «En finir avec le mal de dos» («No more back pain»)**

*«The Climsom Flex heated cervical traction device helps relieve neck pain by encouraging the return to a physiological cervical lordosis.  
It decompresses the vertebrae, relieves the cervical discs, soothes the nervous system and encourages the body's proprioception learning of good postures.  
In synergy, the heating function deeply relaxes the tense muscles in the neck, removing pain, tension and reflex contractures.*

*I highly recommend this effective and pleasant solution, which promotes a liberated neck, capable of healthy and aesthetic postures again.»*

## Cleaning:

Isolate the heating strip from the cover then hand wash the cover with mild soap and water OR use cotton lightly soaked in hydrogen peroxide to rub on the cover.  
More information by flashing the QR Code.

## Composition:

Cervical cushion (SN-FC196): PU leather, density: 35D  
Heating cover (SN-FC196-cover): 14% elastane, 86% nylon-6  
Heating pad: graphene

For more advice,  
flash this QR Code:

www.climsom.com