Mobilly

LAMPE DE LUMINOTHÉRAPIE
Light Therapy Lamp

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
Dispositif Médical

INTRODUCTION

Congratulations on your purchase of this DAYVIA light therapy lamp.

Well used, it will allow you to find the light you lack.

Lack of light is manifested in some people by several of these symptoms:

- Irritability and moodiness
- Significant increase in appetite
- Decreased activity and efficiency at work because of a lack of energy
- Sleep disorders, chronic fatigue
- Problems with sexuality (decreased libido)

Numerous international studies show that light therapy helps restore energy, form, ability to concentrate and good mood like that experienced when you sunbathe in the summer.

The light therapy lamp also helps to reset the waking / sleeping rhythm for people who have shifted work schedules (night workers), people who suffer from phase delay (late falling asleep and difficulty waking up in the morning) or phase advance (too early awakening) or people who travel and cross several time zones (jet-lag).

The effect of light depends on its intensity but above all, on its spectral curve (color).

Indeed, some receptors located in our eyes are very sensitive to light comprising certain wavelengths around 480 nanometers (turquoise blue in particular), they have a major role in our biological clock.

A predominantly blue light may be as effective as a much brighter white light.

Your device therefore offers 2 therapeutic modes, each equally effective:

- intense bluish white light
- a predominantly blue light

We invite you to read this user manual carefully to better benefit from the advantages of your DAYVIA lamp.

SECURITY MEASURES

- The device must be installed and put into service in accordance with the information provided in the accompanying documents.
- Do not use the DAYVIA lamp in a humid environment (for example in a bathroom, near a bathtub or shower, outdoors, etc.).
- Keep the power cable away from heated surfaces.
- Use only the Mean Well power supply box supplied with your light therapy lamp.

This power supply box complies with the requirements for medical devices and is specific for your lamp, it should not be used with other equipment.

- The mains plug acts as a power disconnect switch. The device must be used with the mains plug easily accessible.
- This lamp must not be used by children or people with reduced physical or intellectual capacities.
- Any modification of the device is prohibited.
- The device must not be used if the power cable and / or the power supply unit are apparently damaged or defective.
- The lamp needs special precautions, it can be installed in all premises including domestic premises and away from HF sources where the intensity of electromagnetic interference is high.
- Do not use this device next to other devices or stacked with them because this may cause a malfunction. If this use is necessary, this device and other devices should be observed to verify normal operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) from any part of the DAYVIA lamp. Otherwise, the performance of these devices could be impaired.

PRECAUTIONS

The light therapy lamp must be used in an environment normally lit by another light source to avoid too much contrast between the lamp and its environment, the lamp is complementary to normal lighting (electric or natural light).

Unlike daylight or other lamps on the market, this medical light therapy device does not emit UV rays.

Studies have shown that light therapy can be performed safely (for comparison, your light therapy lamp diffuses either an intense bluish white light up to 10,000 lux or a specific predominantly blue light, while on a beach, in summer, the illumination can reach up to 100,000 lux with the addition of harmful rays such as UV which your lamp does not have).

However, before starting your treatment, please consult your doctor if:

- You suffer or have suffered from depression
 - You are very sensitive to light
 - You suffer from eye disease
 - You suffer from mental disorders
 - You are undergoing treatment with psychotropic or photosensitizing drugs.
- If you have any problems (ex headache or eye strain) after using the light therapy lamp please tell your doctor.

It is often sufficient to start treatment with shorter sessions (10 minutes) and gradually increase them each day to the recommended duration.

START-UP / SHUTDOWN

The DAYVIA Mobilly lamp is carefully checked before it is packaged and shipped.

You will find the user manual, the power supply box and your lamp in your packaging.

We advise you to keep your packaging if a return is necessary.

Here are the different steps to start up your DAYVIA light therapy lamp:

1. Insert the female connector of the low voltage cable coming out of the lamp into the male connector located on the power supply box supplied with the lamp.
2. Check that the voltage indicated on the power supply unit matches the local mains voltage, then plug in the device.
3. Turn on the lamp, in an already lit environment, with the ON / OFF switch located on the left side at the top of the lamp. You can then direct the lighting part towards your face for the duration of the light therapy session and then modulate the intensity with the - and + buttons located next to the switch. This first press of the switch triggers an intense bluish

white light that is very effective from a therapeutic point of view.

You can reduce the intensity with the - button, the light will get warmer and warmer for additional lighting. A second press on the ON / OFF triggers a predominantly blue light of lower intensity as effective as the intense bluish white light of mode 1. A third press on the ON / OFF proposes a neutral white light (without blue) for lighting. supplement without therapeutic aim.

4. After the end of use, turn off the device by simply pressing the ON / OFF switch (the shutdown is done after scrolling through the 3 modes, so on the 4th press ON / OFF).

5. Unplug the power cable from the AC outlet.

RECOMMENDATIONS FOR USE

Place the lamp between 30 and 70 cm from the eyes in the field of vision, directing the illuminating part towards the face during the session, no need to look at it continuously with your eyes, continue your activity normally (breakfast, reading, office, computer, makeup...).

For pleasant, glare-free use, always make sure the room is lit (especially the first few days

of use) to minimize the contrast between the lamp's illumination and its surroundings.

It is advisable to start light therapy sessions as soon as you feel a lack of light before the first symptoms of the winter blues appear.

In prevention or treatment of winter blues, it is preferable to practice the sessions in the morning after waking up.

Do one session each morning per 3 week cycle during the dark period. These figures are to be adapted to each person according to the results obtained after a few days of treatment.

Mode 1: Bluish white light option maximum intensity: 30 min at 30 cm / 45 min at 50 - 70 cm


Mode 2: Light with predominantly blue option maximum intensity: 45 min at 30 cm / 60 min at 50 cm

Mode 3: Non-therapeutic white light option:
Non-therapeutic supplemental light

In general, all people wishing to compensate for their lack of light (winter season, work in a confined environment ...) by light therapy should favor an exposure in the morning: We recommend an exposure between 30 minutes and one hour (see above depending on distance) at the maximum intensity level. This morning exposure, ideally just after waking up, is also recommended for people late in the phase (falling asleep too late and difficulty emerging in the morning) and for people preparing 2 to 3 days before a trip to the east (crossing more than 3 time zones), the effects of jet lag will thus be minimized.

Light therapy at the end of the day is reserved for people working at night, people in early phase (going to bed too early and getting up too early in the morning) and people preparing for a trip to the west (crossing more than 3 time zones) to avoid jet lag syndrome.

After your daily light therapy session, you can use your DAYVIA Mobilly lamp as a desk lamp throughout the day.

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Unplug the lamp and the power adapter and allow it to cool before cleaning it.

Never use abrasives or detergents to clean the device: Use a damp cloth, making sure that no liquid enters the device.

Regarding LED technology, there is no maintenance and the lifespan of the LEDs is more than 30,000 hours.

CHARACTERISTICS OF THE DEVICE

Specifications are subject to change without notice.

Power supply characteristics:

- Manufacturer: Mean Well
- Model: GSM25E15-P1J
- Input voltage: 100-240 V
- Input current: 0.7-0.35 A
- Power supply output voltage: 15V
- Output current: 1.66 A
- Power: 25 W
- Frequency: 50-60 Hz
- Insulation: Class II
- Penetration of foreign bodies and water: IP22

Specifications are subject to change without notice.

Power characteristics:

- Manufacturer: Mean Well
- Model: GSM25E15-P1J
- Input voltage: 100-240 V
- Input current: 0.7-0.35 A
- Power supply output voltage: 15V
- Output current: 1.66 A
- Power: 25 W
- Frequency: 50-60 Hz
- Insulation: Class II
- Penetration of foreign bodies and water: IP22

Lamp characteristics:

- Operating temperature: between + 5 ° C and + 40 ° C
- Humidity rate during use: between 15% and 90% (non-condensing)
- Penetration of foreign bodies and water: IP21
- Atmospheric pressure during use: between 700 hPa and 1060 hPa
- Storage and transport temperature: between - 25 ° C and + 70 ° C with relative humidity up to 90% without condensation and water vapor pressure not exceeding 50 hPa
- Overvoltage category II
- Cord length: 3m22

- Product dimension: 34 x 9 x 9 (W x H x D in cm)
- Weight of the device excluding power supply: 690 gr
- Weight of the device with power supply: 910 gr
- Class I medical device

DAYVIA Mobilly has been certified to comply with the requirements of Directive 93/42 / EEC and 2007/47 / EEC applicable to medical devices in Europe:

- Electrical safety compliance
- Compliance in terms of photobiological risks of luminaires: Risk group 0 (no risk)
- Electromagnetic compatibility compliance: Compliance for each standard or emissions and immunity tests specified by standard EN 60601-1-2: 2015. Class B and group 1 (RF emissions CISPR 11). Immunity to electrostatic discharges: Discharges in air 15kV, on contact 8 kV. Immunity to electromagnetic fields radiated at radio frequencies: 10V / m. Immunity to fast electrical transients in bursts: 2 kV. Immunity to shock waves: ± 1 kV differential mode and ± 2 kV common mode. Immunity to conducted disturbances induced by radioelectric fields: 3V from 150kHz to 80MHz, 6V from for other frequency bands.

Immunity to the magnetic field at the network frequency: 3A / m. Immunity to short interruptions and voltage dips: 0% 10ms, 70% 500ms, 0% 5000ms.

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EXPLANATION OF SYMBOLS

Consult the instructions for use provided.

Manufacturer

Catalog reference

Serial identification

This device is class II

Year of manufacture

This device should not be thrown away with household waste when it reaches the end of its life, drop it off at a collection point where it can be recycled. In doing so, you will help protect the environment.

IP21 Degrees of protection provided by enclosures against liquids and particles.

Limit temperatures


Humidity limits

This device is certified to comply with the Medical Devices Directive.

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