**DAYVIA®**  
LUMINOTHERAPIE

# Wellbi

LAMPE DE LUMINOTHÉRAPIE  
*Light Therapy Lamp*

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Dispositif Médical

## INTRODUCTION

Congratulations on your purchase of this DAYVIA light therapy lamp.

Well used, it will allow you to find the light you lack.

Lack of light is manifested in some people by several of these symptoms:

- Irritability and moodiness
- Significant increase in appetite
- Decreased activity and efficiency at work because of a lack of energy
- Sleep disorders, chronic fatigue
- Problems with sexuality (decreased libido)

Numerous international studies show that light therapy helps restore energy, form, ability to concentrate and good mood like that experienced when you sunbathe in the summer.

The light therapy lamp also helps to reset the waking / sleeping rhythm for people who have shifted work schedules (night workers), people who suffer from phase delay (late falling asleep and difficulty waking up in the morning) or phase advance (too early awakening) or people who travel and cross several time zones (jet-lag).

The effect of light depends on its intensity but above all, on its spectral curve (color).

Indeed, some receptors located in our eyes are very sensitive to light comprising certain wavelengths around 480 nanometers (turquoise blue in particular), they have a major role in our biological clock.

A predominantly blue light may be as effective as a much brighter white light.

Your lamp offers 2 therapeutics modes as effective as each other :

White blue enriched

Cyan Blue targeted to 480nm, the most efficient wave length

## SECURITY MEASURES

-The device must be installed and put into service in accordance with the information provided in the accompanying documents.

-Do not use the DAYVIA lamp in a humid environment (for example in a bathroom, near a bathtub or shower, outdoors, etc.).

- Keep the power cable away from heated surfaces.

-Use only the Mean Well power supply box supplied with your light therapy lamp. This power supply box complies with the requirements for medical devices and is specific for your lamp, it should not be used with other equipment.

-The mains plug acts as a power disconnect switch. The device must be used with the mains plug easily accessible.

-This lamp must not be used by children or people with reduced physical or intellectual capacities.

- Any modification of the device is prohibited.

-The device must not be used if the power cable and / or the power supply unit are apparently damaged or defective.

-The lamp needs special precautions, it can be installed in all premises including domestic premises and away from HF sources where the intensity of electromagnetic interference is high.

-Do not use this device next to other devices or stacked with them because this may cause a malfunction. If this use is necessary, this device and other devices should be observed to verify normal operation.

-Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) from any part of the DAYVIA lamp. Otherwise, the performance of these devices could be impaired.

## PRECAUTIONS

The light therapy lamp must be used in an environment normally lit by another light source to avoid too much contrast between the lamp and its environment, the lamp is complementary to normal lighting (electric or natural light).

Unlike daylight or other lamps on the market, this medical light therapy device does not emit UV rays.

Studies have shown that light therapy can be performed safely (for comparison, your light therapy lamp diffuses either an intense bluish white light up to 10,000 lux or a specific predominantly blue light, while on a beach, in summer, the illumination can reach up to 100,000 lux with the addition of harmful rays such as UV which your lamp does not have).

However, before starting your treatment, please consult your doctor if:

- You suffer or have suffered from depression
- You are very sensitive to light
- You suffer from eye disease
- You suffer from mental disorders
- You are undergoing treatment with psychotropic or photosensitizing drugs

If you have any problems (ex headache or eye strain) after using the light therapy lamp please tell your doctor. It is often sufficient to start treatment with shorter sessions (10 minutes) and gradually increase them each day to the recommended duration.

## START-UP / SHUTDOWN

The DAYVIA Wellbi lamp is carefully checked before it is packaged and shipped.

You will find the user manual, the power supply box and your lamp in your packaging.

We advise you to keep your packaging if a return is necessary.

Insert the removable bracket in the two holes provided on the bottom of the back of the lamp.


You also have the possibility of opting for a wall hanging: use the strap on the top of the lamp and hang it on a screw on the wall.

Here are the different steps to start up your DAYVIA light therapy lamp:


1. Insert the male connector of the power cord supplied with the lamp into the female connector of the low voltage cord coming out of the lamp.
2. Check that the voltage indicated on the power supply unit matches the local mains voltage, then plug in the device.

3. The lamp controls :

4. Turn on the lamp, in an already lit environment, with the on / off switch located at the bottom right of the lamp (middle button above : 1). This first press of the switch triggers an intense bluish white light (default mode) that is very effective from a therapeutic point of view.

5. The time selection button (3) is used to set the time of your light therapy session, after which the lamp switches off automatically.

-Once the lamp is on the White + Blue light Mode, the standard time is 30 min with the high intensity.

-A press on 3 must pass the lamp into medium intensity for 45 min

-Another press on 3 must pass the lamp into low intensity for 60 min.

6. Then a press on the mode selection button above (2) triggers a predominantly blue light of lower intensity as effective as the intense white bluish light.

-Once the lamp is on the Blue light mode, the standard time is 30 min with the high intensity.

-A press on 3 must pass the lamp into medium intensity for 45 min

-Another press on 3 must pass the lamp into low intensity for 60 min.

7. Another press on the mode selection button (2) offers a neutral white light (without blue) for additional lighting without therapeutic aim. The standard is the high intensity.

- A press on 3 must pass the lamp into medium intensity.

-Another press on 3 must pass the lamp into low intensity.

There is not Timer on this light mode.

-After that, another press on 2 bring us to the bluish white light (default mode) and the cycle starts all over again.

8. For the therapeutical mode, the time decreases, and the indicator lights go out progressively, once the programmed time has been used up, the lamp switches off automatically.

You can turn off the device by simply pressing the on / off switch (1).

## RECOMMENDATIONS FOR USE

Place the lamp between 30 and 70 cm from the eyes in the field of vision (approximately 30 ° angle in front of your face), no need to look at it continuously with your eyes, continue your activity normally (breakfast, reading, office work , computer, makeup...).

For pleasant, glare-free use, always make sure the room is lit (especially the first few days of use) to minimize the contrast between the lamp's illumination and its surroundings.

It is advisable to start light therapy sessions as soon as you feel a lack of light before the first symptoms of the winter blues appear.

In prevention or treatment of winter blues, it is preferable to practice the sessions in the morning after waking up.

Do one session each morning per 3 week cycle during the dark period. These figures are to be adapted to each person according to the results obtained after a few days of treatment.

In general, all people wishing to compensate for their lack of light (winter season, work in a confined environment ...) by light therapy should favor an exposure in the morning: We recommend an exposure time of between 30 and 60 minutes.. This morning exposure, ideally just after waking up, is also recommended for people late in the phase (falling asleep too late and difficulty emerging in the morning) and for people preparing 2 to 3 days before a trip to the east (crossing more than 3 time zones), the effects of jet lag will thus be minimized.

Light therapy at the end of the day is reserved for people working at night, people in early phase (going to bed too early and getting up too early in the morning) and people preparing for a trip to the west (crossing more than 3 time zones ) to avoid jet lag syndrome.

After your daily light therapy session, you can use your DAYVIA Wellbi lamp as a desk lamp throughout the day by directing the lighting part towards the worktop and reducing the intensity of the white light in mode 1 (first press on on / off) or by opting for mode 3 (3 presses on / off) according to your wishes.

## DEVICE MAINTENANCE

Unplug the lamp and the power adapter and allow it to cool before cleaning it.

Never use abrasives or detergents to clean the device: Use a damp cloth, making sure that no liquid enters the device.

Regarding LED technology, there is no maintenance and the lifespan of the LEDs is more than 30,000 hours.

## CHARACTERISTICS OF THE DEVICE

Specifications are subject to change without notice.

Power supply characteristics:

- Manufacturer: Mean Well
- Model: GSM25E24-P1J
- Input voltage: 100-240 V
- Input current: 0.7-0.35 A
- Power supply output voltage: 24V
- Output current: 1.04 A
- Power: 25 W
- Frequency: 50-60 Hz
- Insulation: Class II
- Penetration of foreign bodies and water: IP22

Lamp characteristics:

- Operating temperature: between + 10 ° C and + 40 ° C
- Humidity rate during use: between 20% and 90% (non-condensing)
- Atmospheric pressure during use: between 700 hPa and 1060 hPa
- Storage and transport temperature: between - 30 ° C and + 75 ° C with relative humidity up to 90% without condensation and water vapor pressure not exceeding 50 hPa
- Overvoltage category II
- Cord length: 2 m
- Product dimension: 30 x 46 x 1.7 (W x H x D in cm)
- Weight of the device excluding power supply: 2.2kg
- Weight of the device with power supply: 2.4 kg
- Class I medical device

DAYVIA WELLBI has been auto-certified to comply with all requirements of the European Medical Device Regulation (EU) 2017/745 and Directive 2007/47/EEC applicable to medical devices in Europe:

- Electrical safety compliance
- Compliance in terms of photobiological risks of luminaires: Risk group 0 (no risk)
- Electromagnetic compatibility compliance: Compliance for each standard or emissions and immunity tests specified by standard EN 60601-1-2: 2015. Class B and group 1 (RF emissions CISPR 11). Immunity to electrostatic discharges: Discharges in air 15kV, on contact 8 kV. Immunity to electromagnetic fields radiated at radio frequencies: 10V / m. Immunity to fast electrical transients in bursts: 2 kV. Immunity to shock waves: ± 1 kV differential mode and ± 2 kV common mode.

Immunity to conducted disturbances induced by radioelectric fields: 3V from 150kHz to 80MHz, 6V from for other frequency bands. Immunity to the magnetic field at the network frequency: 3A / m. Immunity to short cuts and voltage dips: 0% 10ms, 70% 500ms, 0% 5000ms.

## GUARANTEE AND ASSISTANCE

If you need further assistance or information, visit the DAY-VIA at <http://www.dayvia.com/>

## EXPLANATION OF SYMBOLS

Consult the instructions for use provided.

Manufacturer

Catalog reference

Serial identification

This device is class II

Year of manufacture

This device should not be thrown away with household waste when it reaches the end of its life, drop it off at a collection point where it can be recycled. In doing so, you will help protect the environment.

**IP22**

Degrees of protection provided by enclosures against liquids and particles.

Limit temperatures

Humidity limits

Medical Devices

This device is certified to comply with the Medical Devices Directive.

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DAYVIA SAS