

INNOLUX MESA MEGA LED

Product number: 703280L
(6420611984805)



Certified medical device

EN ISO 13485 & CE 2460

SIZE: k 635, l 405, s 110 / 5,3 kg

THERAPY DISTANCE: 140 cm (2500 lx),
70 cm (10 000 lx)

POWER: 85W LED

LAMP BASE: Fixed LED, 4000K

MATERIAL: Internal painted steel

CORD: acrylic 300 cm white cord

Stepless dimmer

INNOLUX MESA MEGA LED is one of the most powerful bright light devices on the market. The wide light surface combined with high-quality LED components provides a smooth and flicker-free light therapy. The stepless dimmer button enables versatile usage of the lamp. Mesa Mega LED is more energy efficient and has better lighting features than its previous fluorescent light model.

INNOLUX