

INNOLUX SUPERNOVA LED

Product number: 703480L
(6420611984812)



Certified medical device

EN ISO 13485 & CE 2460

SIZE: k 635, l 405, s 110

WEIGHT: 5,5 kg

THERAPY DISTANCE: 200cm (2500lx),
85 cm (10 000 lx)

POWER: 160W LED

LAMP BASE: Fixed LED, 4000K

MATERIAL: painted metal, polycarbonate

CORD: 300 cm white cord

Stepless dimmer

INNOLUX SUPERNOVA LED is made from carefully selected domestic components and suitable for demanding bright light treatment. The truly effective therapy distance ensures the benefits of bright light treatment up to 200 cm away from the lamp. The stepless dimmer enables versatile usage of light. Flicker-free and more even distribution of light guarantees high quality bright light therapy.

INNOLUX